


Cell viability assays with Mithras LB 940

Emmanuelle Estace¹, Veronika Keravec¹, Bettina Felletschin²

Introduction

Actinomycin D (picture 1) is a chemical compound with the formula $C_{62}H_{86}N_{12}O_{16}$ and molecular weight of 1255.41 g/mol. It is used as an anti-mitotic and an intercalating agent of DNA. It inhibits the synthesis of mRNA by intercalation into the DNA strand and thereby preventing the polymerase activity.

In this study we tested the effects of two concentrations of Actinomycin D on CHO cells. The cellular viability of those cells was assessed by different cell viability agents/assay kits like Alamar Blue (Invitrogen), UptiBlue (Interchim) and Resazurin (Interchim).



Picture 1 : Actinomycin D

The assays are easy-to-use as no cell extraction or lysis steps are required in contrast to other commonly used cell proliferation assays. All indicators are non-toxic and water soluble thus eliminating washing/fixing steps.

The systems are based on a dye (Resazurin) which is highly sensitive to the redox potential of cell growth. The dye changes by responding to the chemical reduction of growth medium from the oxidized (non-fluorescent, blue) to the reduced (fluorescent, red) form.

¹ INSERM Unité 644 - Medicine and Pharmacy Institute – University of ROUEN, 76183 ROUEN Cedex France

² BERTHOLD TECHNOLOGIES, Bad Wildbad, Germany, Bettina.Felletschin@Berthold.com

The Mithras LB 940 is a multimode plate reader with an unique optical design (DOPS – Dedicated Optical Path System) to ensure optimized performance for the detection technologies implied. These are

- luminescence
- BRET/BRET²
- fluorescence
- UV/VIS absorbance
- fluorescence polarization
- AlphaScreen[®]
- TRF
- HTRF[®]




Picture 2: Mithras LB 940 multimode reader

In addition accessory options, e.g. reagent injectors, temperature control and cooled PMT detection units are available.


Instrument Settings

Following settings in MikroWin 2000 were used:




Alternatively a filter combination of 530/10nm for excitation and 600/10nm for emission can be used as well.

Triplicates were selected in the Plate Partition part:




Under Options/Results both the "Reader" and "Sample" matrix were added to the "Results" matrix to get the sample (name) information, single well data and average of the three wells displayed into a single screen:



The screenshot shows the MikroWin software interface with a data table and an 'Options Result' dialog box. The data table has columns for wells 1-9 and rows A-H. The 'Options Result' dialog box is open, showing a 'Matrix' dropdown set to 'Results' and a 'Combination' table with two entries: 'Reader No. 1' and 'Sample No. 2'. There are 'Add Matrix' and 'Delete Matrix' buttons in the dialog.

#	1	2	3	4	5	6	7	8	9
A	0	0	0 Control	0	0	0 Act D 0.5 µM	0	0	0 Act D 2 µM
B	0						0	0	0 Act D 2 µM
C	0						0	0	0 Act D 2 µM
D	0						0	0	0 Act D 2 µM
E	0						0	0	0 Act D 2 µM
F	0						0	0	0 Act D 2 µM
G	0	0	0 Control	0	0	0 Act D 0.5 µM	0	0	0 Act D 2 µM
H	0	0	0 Control	0	0	0 Act D 0.5 µM	0	0	0 Act D 2 µM


The print setup was set as the following:



Measurement file : upti 120908 3h.dat Valid Assay Measurement date : 12/09/2008
 Template file : UptiBlue Test Berthold.par Measurement time : 13:44:15

Echantillons / Données brutes / Moyenne												
	1	2	3	4	5	6	7	8	9	10	11	12
A	41160,0	38220,0	38310,0	39230,0	40660,0	38380,0	37440,0	38826,7	39940,0	37210,0	37270,0	38140,0
B	37960,0	36620,0	36610,0	37053,3	41700,0	38380,0	38230,0	39436,7	39910,0	39290,0	37350,0	38650,0
C	712680	690880	682510	695357	479500	448850	447000	458450	509710	482670	490530	494303
D	1,04E6	995900	997280	1,01E6	425400	393080	398810	405763	414950	379490	382760	392400
E	799860	736130	739650	758547	504740	469000	471390	481710	474460	450620	445920	457000
F	743230	749960	743980	745723	433490	410340	408300	417377	425530	401140	400410	409027
G	702100	704980	703440	703507	512470	475570	478840	488960	402190	386120	367700	378670
H	726930	727620	731460	728670	413690	384800	383600	394030	420880	387340	391810	400010

General Data
 Measurement EndPoint : Time of Measurement : 12/09/2008 / 13:44:15, State of Measurement : Valid Measurement
 Parameter Parameter loaded during measurement time : c:\program files\mikrowin 2000\paramathras\fluoresc
 Files Data file : upti 120908 3h - 12/09/2008 / 12:44:48
 Template file : UptiBlue Test Berthold - 05/09/2008 / 10:45:20
 Reader BertholdTech Mitthras, Driver Version: 1.05 , (1.0.5.0), S/W: 33-6002, Embedded Version: 1.17
 Plate Type: 8x12 plate
 Name Fluor. Label
 Counting Time [s] 0.10 Measurement Mode by Plate
 Excitation Filter 520 Emission Filter 590
 Excitation Aperture Normal
 Lamp Energy 7000 Counter Position Top
 Second Measurement No
 Calculation Calculation Status : Valid Assay, Calculation Time : 16/12/2008 / 09:32:44
 Program MikroWin, Version 4.40; License No. : 2997; Assembly Code : 0066 0001 FFFB FFFB
 Operating System : Windows 2000 Professional; User Name : Toxicologie; Printer Name : HP Laser.



Methods

CHO cells were grown in IMDM medium (Invitrogen) supplemented with 10 % FCS (Fetal Calf Serum, Invitrogen), HT supplement (SIGMA) and penicillin/streptomycin (Invitrogen). The cells were cultivated in 75 cm² flasks (Falcon) in a humid atmosphere incubator at 37 °C under 5 % CO₂ to produce enough cells and then, after trypsination, seeded in 6 wells plates (Falcon) to perform the experiments. The medium was changed every 2 days until testing process.

Actinomycin D was prepared in IMDM medium and used at two concentrations (0.5 and 2 µM) versus control (IMDM medium alone).

Alamar Blue (Invitrogen), UptiBlue (Interchim) or Resazurin (Interchim) was added each 150 µL to the 6 well plates. The fluorescence readings were performed in 96 well plates to avoid toxic effects and less oxygenation. A volume of 100 µL was transferred from the 6 well plates to the 96 well black plates (Nunc) to perform the readings with the Mithras LB 940 (Berthold Technologies).

As indicated in the kit instructions, the incubation time had to be evaluated according to the cell type. The incubation of CHO cells with UptiBlue and Resazurin was done then for 3 hours.

Cell density determination

First, cell density experiments were performed to determine the optimum cell numbers. Four 6 wells plates were seeded with different cell densities: 0, 50 000, 100 000, 300 000 and 500 000 cells per well in 1.5 mL IMDM medium. After 2 days the cell culture medium was discarded, and the cellular monolayer was washed with 1 mL Ca²⁺/Mg²⁺ free PBS solution. Finally, 1.5 mL IMDM fresh medium and 150 µL UptiBlue, Alamar Blue or Resazurin solution (dilution 1/10e) were added to each well according to the manufacturer's instructions. Cells were kept in a humid atmosphere incubator at 37 °C and 5 % CO₂ for 3 hours. At the end of incubation, 100 µL medium was taken in duplicate from each well and transferred into 96 well black plates for reading.

The readings were performed with Mithras LB 940 multimode reader and Mikrowin 2000 software. A fluorescence parameter file with an excitation wavelength at 520 nm and an

emission wavelength at 590 nm and a counting time of 0.1 s per well was used as described before in instrument settings.

Cytotoxicity Assay


The cytotoxicity assay was performed in 6 well plates with 100 000 cells per well with the same medium and volume as used before. After 4 days cell culture, the medium was removed and each well was replaced with fresh medium containing either Actinomycin D at 0.5 μ M and 2 μ M or fresh medium without any Actinomycin D (control). Plates containing all reagents but no cells were used as a negative control to determine the potential of non-cell-related reduction of UptiBlue, Resazurin or Alamar Blue over the experimental period. After 21 hours incubation 150 μ L Alamar Blue, Resazurin or UptiBlue solution was added to each well after medium removal procedure or not.

For medium removal procedure fresh medium (without any Actinomycin D) was added to each well after elimination of the old medium by aspiration. The 150 μ L fluorophore solution was added into this fresh medium. Without medium removal procedure the fluorophore solution was added directly into the old medium.

Then, incubation was continued for 3 hours under the same conditions. At the end of incubation, the same reading protocol and transfer to 96 well plates was performed as used for the cell density determination.


For protein quantification, the medium containing the fluorophore was discarded after the reading. The cells were washed with 2mL of Ca²⁺/Mg²⁺ free PBS solution. Then 0.5mL of NaOH 1N solution was added to dissolve protein and to quantify protein content according to Lowry and al. (1951). The absorbance readings were performed with the Mithras at 650 nm.

Results




Picture 3: Results of cells density determination, relative fluorescence units according to the cell number


Picture 3 shows that the UptiBlue, AlamarBlue or Resazurin reduction was exactly proportional to the cell seeding number. With more than 300 000 cells the response approaches a saturation level. The reason may be either a depletion of fluorescent substrate and/or metabolic stress related to overcrowding.




Picture 4 : Photographs of cell culture (with or without medium removal procedure)



Picture 5 : Comparison of UptiBlue, Alamar Blue and Resazurin reduction without medium removal procedure




Picture 6 : Comparison of UptiBlue, Alamar Blue and Resazurin reduction with medium removal procedure



Picture 7 : Comparison of attached cell protein content

Correlation between fluorophore reduction and attached cell protein content



Without medium removal procedure, neither by using UptiBlue, Alamar Blue or Resazurin differences between Actinomycin D solutions and control could be observed.

However, after medium change (medium removal procedure), a strong decrease of UptiBlue, Alamar Blue or Resazurin reduction in the Actinomycin D treated wells could be observed.

The CHO cells are specifically able to get suspended when they are in stress. In our case this metabolic stress was caused by the presence of Actinomycin D and the cells got suspended in culture medium.

According to the obtained results a part of the cells were still active. Finally, in our experimental conditions, it was necessary to discard medium containing suspended cells and add fresh medium before performing the cytotoxicity assay.

All different reagents used (UptiBlue, AlamarBlue and Resazurin) have shown comparable results.

Moreover, these experiments have shown showed a good correlation between the

chromophore reduction and the seeded cell numbers and the protein content of the attached cells.

Differences between the two Actinomycin D concentrations could be observed by the use of Resazurin as indicating reagent only. Additional measurements are necessary to get more detailed information.

Conclusion

Alamar Blue, UptiBlue and Resazurin are perfectly suited for performing cell viability and cytotoxicity assays in combination with the Berthold Mithras microplate reader.

The medium removal procedure before performing the cytotoxicity assay is essential to avoid influences of suspended cells and old medium.

Material:

- Mithras LB 940 (Berthold Technologies)
- Filters (Berthold Technologies):
Excitation 520/10nm (39802) or Excitation 530/10nm (37996)
Emission 590/20nm (37989) or Emission 600/10nm (40095)
Absorbance filter 650/10nm (37999)
- CHO cells (Chinese Hamster Ovary)
- IMDM medium (Invitrogen)
- FCS (Fetal Calf Serum, Invitrogen),
- HT supplement (SIGMA)
- Penicillin/streptomycin (Invitrogen)
- 75 cm² flask (Falcon)
- Alamar Blue SKU# DAL1025 (Invitrogen)
- UptiBlue UP669412 (Interchim)
- Resazurin (Interchim)
- 96 wells black plates (Nunc)
- 6 wells plates (Falcon)